

# Analysis of THC and CBD Content in Hemp Oil Using GC-FID and Hydrogen Carrier Gas

Yuan Lin, Shimadzu Scientific Instruments

## Abstract:

Hemp oil containing beneficial phytocannabinoids such as cannabidiol (CBD) but very low levels of the psychoactive compound tetrahydrocannabinol (THC) has become increasingly popular. For quality control of the product, the level of CBD and THC should be closely monitored. In this report, hemp oil samples were analyzed by the Shimadzu GC-2030 using H<sub>2</sub> carrier gas and the concentrations of CBD and THC were determined.

## 1. Introduction

Both hemp and marijuana are varieties of the cannabis plant. While marijuana generally possesses high levels of the psychotropic tetrahydrocannabinol,  $\Delta^9$ -THC, and lower levels of the non-psychotropic cannabidiol, CBD, hemp is just the opposite: rich in CBD but low in  $\Delta^9$ -THC. While more research is needed to better understand the chemistry of benefits from CBD, it has been reported to reduce convulsions, inflammation, nausea and anxiety, and has even eradicated tumors in some patients. (Note: the cannabinoid(s) recommended for specific medical conditions have not been approved by the FDA.)

CBD-rich oil has become increasingly popular and is administered via sublingual drops, gel capsules or as a topical ointment. The main source of CBD-rich oil is industrial hemp. CBD oil is derived as concentrate from CO<sub>2</sub> or butane extraction of hemp, sometimes followed by steam distillation or ethanol distillation for purification.

The FDA has issued warning letters to firms that market unapproved new drugs allegedly containing CBD. As part of these actions, the FDA has determined the cannabinoid content of some hemp products and many were found to contain levels of CBD that are very different from the label claim. It is important to note that such products are not approved by the FDA for the diagnosis, cure, mitigation, treatment, or prevention of any disease.

Therefore, as quality control, the level of CBD and THC in any such products should be monitored. Like cannabis, hemp oil may be analyzed easily and effectively for its cannabinoid content. In this report, both THC and CBD levels in two different hemp oil samples were determined using a GC-FID setup.

## 2. Experimental Methods

The cannabinoids standard of CBD,  $\Delta^9$ -THC and CBN (cannabinol) was purchased from Restek (cat no. 34014). Hemp oil samples were obtained from various vendors.

To make calibration standards, original stock of a cannabinoids standard was diluted in methanol to indicated concentrations. The hemp oil samples were diluted as follows: 10 $\mu$ L of oil was first diluted into 400 $\mu$ L isopropanol, mixed thoroughly, then diluted again into 400 $\mu$ L methanol and mixed thoroughly.

### Analytical Conditions

<b>Column</b>	SH-Rxi-5Sil MS, 15m x 0.25mm x 0.25 $\mu$ m (part no. 227-36036-01)
<b>GC oven temp</b>	200°C ramp at 15°C/min to 300°C, hold 5min
<b>SPL</b>	250°C; Hydrogen carrier gas, constant inlet pressure at 48.3kPa, split ratio = 10
<b>FID</b>	300°C; Hydrogen flow rate 32mL/min; Air flow rate 200mL/min; Nitrogen makeup gas flow rate 24mL/min.

Data were acquired and analyzed using Lab Solutions software.




GC-2030 FID with AOC-20i Autoinjector


### 3. Results and Discussion

A six-point calibration curve was generated using the cannabinoid standards of 0.8ppm, 4ppm, 20ppm, 100ppm, 250ppm and 500ppm each of CBD,  $\Delta^9$ -THC and CBN in methanol and fitted to linear regression (fig.1). The correlation coefficient ( $r^2$ ) values of all three compounds are  $>0.999$  over 6 standard levels.

#### 1. CBD (cannabidiol)



#### 2. THC ( $\Delta^9$ -tetrahydrocannabinol)



#### 3. CBN (cannabinol)




Figure 1: Standard curves for cannabinoids




Figure 2: Example chromatogram of hemp oil sample 1. Note that CBN was not detected in the sample.

	CBD (total)		THC (total)	
	conc (%)	%RSD	conc (%)	%RSD
Sample 1	2.182	0.893	0.047	1.596
Sample 2	0.062	0.555	<i>Not detected</i>	--

Table 1: Measured concentrations of CBD and THC in two different hemp oil samples.

Each oil sample was diluted as described in Experimental Methods and run as triplicates. CBN was not detected in either of these samples. To calculate the % level of cannabinoids in the samples, the reported ppm concentration of the sample from LabSolutions software was multiplied by 81 (dilution factor) then divided by 10,000 (10,000ppm = 1%).

## 4. Conclusions

In this report, cannabinoid CBD and THC levels in two hemp oil samples were analyzed by the Shimadzu GC-2030 with FID and AOC-20i autoinjector. Calibration was carried out from 0.8ppm to 500ppm for each cannabinoid with excellent linearity. No carryover was observed with this setup and method. H<sub>2</sub> was used as the carrier gas for efficient and low cost analysis. Furthermore, with an integrated hydrogen sensor option for GC-2030, the instrument can be used with hydrogen carrier gas routinely with minimal safety concerns for the operator.

As shown in the table of results above, sample 1 is shown to contain a high level of CBD (2.18%), while sample 2 has very little CBD in comparison (0.06%). Clearly not all products are created equal. And this signifies the importance of assaying CBD level in these claimed CBD-rich products. It should be noted, though, that the CBD amount reported in the sample is the summation of CBD and CBDA (Cannabidiolic acid), or total CBD, as the CBDA present would have been converted to CBD in the hot injection port of the GC and added to the CBD peak.

Also shown by the results above, both hemp oil samples contain very low levels of the psychoactive compound THC. The low level of THC (< 0.3%) is what distinguishes hemp from marijuana; therefore, the THC level is another important parameter to monitor in these products. Similar to CBD, one consideration is that the acid form of THC known as THCA will be converted to THC in the hot injection port of the GC; as a result, the THC level reported above is the sum of THCA and Δ9-THC, or total THC.

Furthermore, this method can be adapted easily to assay the THC level in medical marijuana for potency. However, as mentioned above, the level determined will be the sum of THCA and Δ9-THC. If individual compound concentration is desired, an HPLC method should be employed instead. Information regarding Shimadzu HPLC cannabis analyzer can be found on the Shimadzu website [ssi.shimadzu.com](http://ssi.shimadzu.com).

First Edition: December 2018



**Shimadzu Corporation**  
[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

**Shimadzu Scientific Instruments**  
7102 Riverwood Drive, Columbia, MD 21046, USA  
Phone: 800-477-1227/410-381-1227  
[www.ssi.shimadzu.com](http://www.ssi.shimadzu.com)

**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purposes without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks or trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "@". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.