

Technical Report

Oxygen heterocyclic compounds profile evaluation in cold-pressed *Citrus* essential oils using supercritical fluid chromatography coupled to photodiode array detector

Antonella Satira¹, Marina Russo², Paola Dugo^{1,2}, Luigi Mondello^{1,2}

Abstract:

A large number of analytical techniques to characterize the non-volatile composition of *Citrus* essential oils (EOs) were introduced in the literature. However, the reports of *Citrus* EOs analysis by supercritical fluid chromatography (SFC) with UV or mass spectrometry detectors for the analysis of non-volatile polar compounds are still few. The aim of the present scientific research is to describe a new analytical approach using SFC with a green mobile phase for the analysis of such compounds in cold-pressed *Citrus* EOs, to confirm the potential use of SFC-UV for oil classification in the context of quality control of raw materials in cosmetics. Gradient conditions are determined to achieve a satisfactory separation in 10 minutes.

Keywords: Oxygen heterocyclic compounds, *Citrus* essential oils, supercritical fluid chromatography, furocoumarins

1. Introduction


Citrus EOs are cold-extracted from the peel of *Citrus* fruits using different mechanical systems and are composed by a volatile fraction (85–99% of the *Citrus* oils) and a non-volatile fraction (1–15%). Coumarins, psoralens, and polymethoxyflavones are a class of oxygen heterocyclic compounds present in the non-volatile fraction of cold-pressed *Citrus* EOs. It is well known that psoralens exhibit strong photoactivity in combination with UVA radiation [1] and for this reason, there are opinions and regulations on the use of essential oils in cosmetic products [2-4]. The study of the non-volatile fraction of cold-pressed *Citrus* EOs is a valid tool to exploit their quality and authenticity, and it also represents an important means for product control in the cosmetic industry. RP-HPLC represent, from 2000 up to now, the most used technique for the analysis of oxygen heterocyclic compounds with optimal results in terms of accuracy, repeatability, limit of detection, and quantification [5]. In perspective of green chemistry, supercritical fluid chromatography (SFC) was not sufficiently exploited, more investigation using this analytical technique are needed.

The present research is focused on the development of a rapid analytical method with low environmental impact using SFC coupled to a photodiode detector (PDA), to analyze oxygen heterocyclic compounds in cold-pressed *Citrus* EOs. The separation of all compounds was achieved in less than 10 min, using a HILIC stationary phase and a mobile phase including CO₂ and ethanol with 0.6% of H₂O and 0.07% of trifluoroacetic acid (TFA).



Supercritical Fluid Chromatograph Nexera™ UC

AGREE Analytical GREENness Metric Approach



1. Sample treatment and preparation
2. Sample size and number of samples
3. In Situ measurements or in laboratory
4. Analytical steps
5. Automation and miniaturization of analytical method
6. Use of derivatization reagents
7. Volume of analytical waste generated
8. Number of analytes determined
9. Energy consumed
10. Use of reagents obtained from renewable sources
11. Use of toxic reagents
12. Safety of the operator

Figure 2: AGREE graphic of the RP-HPLC-PDA and SFC-PDA analytical methods adopted for the oxygen heterocyclic compounds analysis in cold-pressed *Citrus* EOs.

Table 2: Oxygen heterocyclic compounds identified in cold-pressed *Citrus* EOs investigated using SFC-PDA method.

#	Compounds	Lime	Lemon	Bergamot	Grapefruit	Sweet orange	Bitter orange	Red mandarin
1	Coumarin							
2	Isoimperatorin	x	x					
3	Herniarin	x		x				
4	5-Isopentenyl-7-methoxy-coumarin		x					
5	Citropten	x	x	x	x			
6	Bergamottin	x	x	x				
7	Cnidicin		x					
8	Cnidilin							
9	Psoralen							
10	5-Geranyloxy-7-methoxy-coumarin	x	x	x				
11	Bergapten	x		x	x		x	
12	Aurapten				x			
13	5-(isopenten-2'-eniloxy)-8-(2',3'-epoxy)-isopentenyl-7-methoxy-psoralen	x	x					
14	5-Geranyloxy-7-methoxy-psoralen							
15	Osthol				x			
16	8-Methoxypsoralen						x	
17	Isopimpinellin	x						
18	Heraclenin							
19	Imperatorin	x	x					
20	Phellopterin		x					
21	Epoxybergamottin				x			
22	Epoxyaurapten				x			
23	8-Geranyloxy-psoralen	x	x					
24	Isomeranzin				x		x	
25	Meranzin				x		x	
26	Byakangelicol	x	x					
27	Oxypeucedanin	x	x					
28	Tangeretin				x	x	x	x
29	Byakangelicin	x	x					
30	Epoxybergamottin hydrate			x			x	
31	Oxypeucedanin hydrate	x	x					
32	Heptamethoxyflavone			x		x	x	x
33	Tetra-O-methyl-scutellarein			x		x		x
34	Nobiletin			x		x		x
35	Hexamethoxyflavone						x	
36	Meranzin hydrate			x			x	
37	Sinensetin			x		x		x

4. Conclusion

The separation of oxygen heterocyclic compounds in cold-pressed *Citrus* EOs is a challenge due to the wide variety of compounds and due to small structural differences between the compounds. The results obtained show that SFC-UV is a perfectly suited method to investigate the essential oil composition, because of the great number of compounds separated in a reduced analysis time (around 10 min), and with a very short time for re-equilibration of the system at the end of the gradient analysis.

References

- 1) European Parliament, Official Journal of the European Union, L 342 (2009) 59, 22.12.2009. Regulation (EC) No. 1223/2009 of the European Parliament and of the Council of 30 November 2009 on cosmetic products (recast). <https://eur-lex.europa.eu>.
- 2) International Fragrance Association (IFRA). IFRA Standard. 48th Amendment – Citrus oils and other furo-coumarins containing essential oils. 2015. <https://ifrafragrance.org>.
- 3) International Fragrance Association (IFRA). IFRA Analytical method – quantitative determination of furo-coumarins by HPLC-DAD. 2013. <https://ifrafragrance.org>.
- 4) “European Parliament, Official Journal of the European Union, L 354 (2008) 34, 31.12.2008. Regulation (EC) No 1334/2008 of the European Parliament and of the Council of 16 December 2008 on flavourings and certain food ingredients with flavouring properties for use in and on foods and amending Council Regulation (EEC) No. 1601/91, Regulations (EC) No. 2232/96 and (EC) No. 110/2008 and Directive 2000/13/EC. <https://eur-lex.europa.eu>.
- 5) Marina Russo, Francesca Rigano, Adriana Arigò, Paola Dugo and Luigi Mondello (2021). Coumarins, Psoralens and Polymethoxyflavones in Cold-pressed Citrus Essential Oils: a Review. *Journal of Essential Oil Research*. DOI: 10.1080/10412905.2020.1857855
- 6) Francisco Pena-Pereira, Wojciech Wojnowski, and Marek Tobiszewski (2020). AGREE-Analytical GREENness Metric Approach and Software. *Analytical Chemistry*. DOI: 10.1021/acs.analchem.0c01887
- 7) Marina Russo, Ivana Bonaccorsi, Rosaria Costa, Alessandra Trozzi, Paola Dugo and Luigi Mondello (2015). Reduced time HPLC analyses for fast quality control of Citrus essential oils. *Journal of Essential Oil Research*. DOI: 10.1080/10412905.2015.1027419



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.