

Comprehensive Two-Dimensional HPLC and Informative Data Processing for Pharmaceuticals and Lipids

HPLC 2015 PSB-MULTI-06

Yoshiyuki WATABE, Tetsuo IIDA, Daisuke NAKAYAMA,
Kanya TSUJII, Saki UEDA, Tadayuki YAMAGUCHI,
Junichi MASUDA, Yoshihiro HAYAKAWA
Shimadzu Corporation, Kyoto, Japan

Introduction

Kakkonto, traditional Chinese drug, is believed to be effective against cold, headache, stiff shoulders, etc. It consists of natural plants and contains many compounds such as ephedrine, glycyrrhizic acid, and cinnamic acid. Nexera-e, comprehensive two-dimensional liquid chromatograph, is useful to separate such a complex matrix. In general, different separation modes are selected for 1st and 2nd dimensions for comprehensive two-dimensional chromatography to obtain high resolution that can not be achieved by each individual separation mode. We tried comprehensive separation using a semi-micro reversed-phase column and neutral pH mobile phase for 1st dimension and an ultra-high speed reversed phase column and acidic pH mobile phase for 2nd dimension. Herbal medicine contains many relatively polar compounds, so pH is an important parameter to modify separation selectivity. We performed differential analysis for commercial Kakkonto products and

determined glycyrrhizic acid contained with photodiode array detector (PDA).

In lipidomics, phospholipids are the attractive targets of LC/MS analysis since lipids are important and essential components of biological membranes. However, conventional HPLC system by a single separation mode performs poorly on biological lipid sample, because it contains various kinds of lipids with common moieties that govern their behavior on column. Nexera-e was capable of characterizing phospholipids both quantitatively and qualitatively when coupled with triple quadrupole and iontrap-TOF type of mass spectrometer (MS/MS and IT-TOF) respectively. Reliable identification of lipid species was performed by acquiring m/z values of related parent and fragment ions at high accuracy with the IT-TOF and matching the data to commercially available data-base.

Experimental




Fig.1 Flow line diagram of "Nexera-e"

Comprehensive Two-Dimensional HPLC and Informative Data Processing for Pharmaceuticals and Lipids

Kakkonto

| | |
|----------------|---|
| 1D Column | : Shim-pack XR ODSII (100 mm L. x 1.5 mm I.D., 2.2 µm) |
| Mobile Phase | : A; 10 mmol/L (sodium) phosphate buffer pH = 6.8 B; acetonitrile |
| Flow Rate | : 0.05 mL/min |
| Time Program | : B Conc. 5% (0 min) → 30% (70 min) → 90% (80 min) → 90% (90 min) → 50% (90.1 min) → STOP (110 min) |
| Column Temp. | : 40 °C |
| Injection vol. | : 2 µL |
| Loop vol. | : 50 µL (Modulation time : 60 sec) |
| 2D Column | : Kinetex (50 mm L. x 3 mm I.D., 2.6 µm) |
| Mobile Phase | : A; 10 mmol/L (sodium) phosphate buffer pH= 2.6 B; acetonitrile |
| Flow Rate | : 2 mL/min |
| Time Program | : Without Auto-gradient; B Conc. 5% (0 min) → 60% (0.75 min) → 5%(0.76 min) → STOP (1 min) With Auto-gradient; Initial.B Conc. 5% (0 min) → 40% (0.75 min) → 5%(0.76 min) → STOP (1 min) Final.B Conc. 20% (0 min) → 60% (0.75 min) → 20%(0.76 min) → STOP (1 min) The initial and final B conc. has been changed by a stepwise method |
| Detector | : SPD-M30A Photo diode array detector (PDA, standard cell 1µL, wavelength = 254 nm) |

Glycerophospholipids

| | |
|----------------|--|
| 1D Colum | : Nucleosil SIL (150 mm L. X 1.0 mm I.D., 3 µm) |
| Mobile Phase | : A; isooctane/acetone/ethyl acetate/acetic acid = 40/20/20/0.03 (v/v/v/v) B; isooctane/2-propanol/water/acetic acid/28% ammonia aq.sol.= 40/51/9/0.03/0.03 (v/v/v/v/v) |
| Flow Rate | : 0.02 mL/min |
| Time Program | : B Conc. 30% (0 min) → 40% (25 min) → 100% (40 min) → 100% (55 min) → 30% (55.1 min) → STOP (70 min) |
| Column Temp. | : 40 °C |
| Injection vol. | : 5 µL |
| Loop vol. | : 20 µL |
| 2D Column | : Phenomenex Kinetex C18 (50 mm L. X 4.6 mm I.D., 2.6 µm) |
| Mobile Phase | : A; methanol/water/acetic acid/28% ammonia aq.sol. = 90/10/0.05/0.05 (v/v/v/v) B; 2-propanol/acetic acid/28% ammonium hydroxide = 100/0.05/0.05 (V/V/V) |
| Flow Rate | : 3.5 mL/min (50% split for MS) |
| Time Program | : B Conc. 10% (0 min) → 50% (0.75 min) → 10%(0.76 min) → STOP (1 min) The initial B conc. has been changed by a stepwise method |
| Detector | : Shimadzu LCMS-8050 (MS/MS, ESI positive, MRM mode) Shimadzu LCMS-IT-TOF (IT-TOF, ESI positive/negative scan, m/z :700-800) |

Results

Kakkonto (Chinese crude drug): RP×RP-PDA

Effect of auto gradient function and determination of Glycyrrhizic acid




Fig.2 Comprehensive-2D separation of commercial Kakkonto product with/without "Auto gradient program" function, linearity of 50-1000 mg/L and repeatability of 5 replicated analyses of real sample

Differential analysis between two kakkonto products




Fig.3 Comparison of comprehensive contour plots of two Kakkonto products

Fig.4 Differential analysis between A and B products

Comprehensive Two-Dimensional HPLC and Informative Data Processing for Pharmaceuticals and Lipids

Conclusion

1. Novel comprehensive two-dimensional LC “Nexera-e” was successfully applied to quantitative, qualitative and differential analyses of complex samples.
2. PDA, MS/MS and IT-TOF detection system could go with comprehensive two-dimensional separation with satisfactory sensitivity and ruggedness.
3. Combination use of different separation modes afforded sufficiently orthogonal separation of Chinese crude drug and phospholipids.
4. Excellent repeatability, linearity of blob area as well as repeatability of total retention time by “Nexera-e”-MS/MS or PDA setup were confirmed.
5. Reliable identification of target lipid in 2D contour plot was achieved by data-base matching analysis based on high resolution IT-TOF data.