

The study of aging processes using *in vitro* glycation experiments applied to low density lipoproteins using MALDI-QIT-TOF-tandem Mass Spectrometry

IMSC 2012 Pwe-170

Omar Belgacem¹, Matthew Openshaw¹, Grazyna Sobal², Katharina Pock³ and Gerald Stübiger⁴

¹Shimadzu, Kratos Analytical Ltd, Manchester, UK;

²Dept. of Nuclear Medicine, Medical University Vienna, Austria.

³Octapharma Pharmazeutika Produktionsges.m.b.H, Vienna, Austria.

⁴Center of Biomolecular Medicine and Pharmacology, Medical University Vienna, Austria.

The study of aging processes using *in vitro* glycation experiments applied to low density lipoproteins using MALDI-QIT-TOF-tandem Mass Spectrometry

Introduction

Due to the steadily increased human lifespan during the last centuries our society is increasingly faced with all sorts of health problems related to “aging processes” that are now causing a huge stress on the healthcare system of any modern country. Today glycation of proteins also known as “Maillard reactions” is increasingly recognized to be related to aging processes involved not only in the development of human disease but also in altering the quality of food and pharmaceutical formulations [1, 2]. In this context, the modification of low-density lipoprotein (LDL) induced by elevated glucose levels in combination with reactive oxygen species (ROS) is considered a major risk factor of

cardiovascular disease (CVD) and atherosclerosis [3, 4]. In order to track the range of modifications at the molecular level we have developed a methodology to specifically detect glycation sites on various standard proteins in order to be able to look at real samples such as human serum albumin (HSA), immunoglobulins (e.g. IgG) and lipoproteins (e.g. LDL) *in vivo* [5]. Our lipoproteomics approach using (LC-)MALDI-QIT-TOF-MS/MS combined with an intelligent software tool (PTM Finder™, Shimadzu, UK) represents a unique platform for studying lipids and modification of amino acids (AA) in proteins on the molecular level.




Fig. 1 Workflow of the MALDI-QIT-TOF-MS/MS based lipoproteomics approach

The study of aging processes using *in vitro* glycation experiments applied to low density lipoproteins using MALDI-QIT-TOF-tandem Mass Spectrometry

Experimental

Lipoproteins were isolated from normal human plasma by preparative ultracentrifugation and modified by diverse ROS (e.g. Cu²⁺, HOCl, etc) and lipid peroxidation products (e.g. MDA, HNE) in the presence or absence of sugars (e.g. glucose). PLs were extracted and directly analyzed by MALDI-MS. The apolipoproteins were digested with endoproteinases (e.g. trypsin), peptides were separated by HPLC and sequence mapping was performed based on MALDI-MS/MS. Mass spectra were recorded by MALDI-TOF-RTOF- and MALDI-QIT-TOF-MS/MS using

AXIMA-Performance and AXIMA-Resonance (Shimadzu, UK), respectively. Lipids were analyzed using THAP in 80/20; MeOH/100mM DAHC (v/v) and for peptides CHCA (5 mg/ml) in 50/50 ; ACN/0.1% TFA (v/v) containing 10mM ADHP was used [6, 7]. Data were processed and searched for specific MS/MS signatures based on a parameterizable database and software module (PTM Finder, Shimadzu). Lipids are identified and analysed using prototype software tools (LipIdent) developed in cooperation with external partners. An overview of our workflow is provided in Fig. 1.

Results

Phospholipids (PLs) were recorded in positive ion MALDI-TOF-MS mode in the *m/z* range 400-900 and MALDI lipid profiles of the different *in vitro* modified LDL particles were established (Fig. 2). A decrease of PC containing sn-2 PUFA (i.e. 18:2, 20:4, 22:6) between 19-75%, and increase of LPC16:0/18:0 between 42-350% compared to nLDL was observed in oxLDL, gLDL, and goxLDL, respectively. This indicates for strong lipid degradation and enzymatic hydrolysis (e.g. Lp-PLA2 activity) during the LDL modification process initiated by ROS and/or glucose. MALDI-TOF-MS/MS spectra of peptides were

acquired and large peptide MS/MS dataset (more than 400 MS/MS spectra per LDL type) were searched for characteristic neutral losses linked to sugar moieties attached to the amino acid backbone (-162, -120 and -36 Da) using the PTM Finder tool (Shimadzu, UK) (Fig. 3 and 4). Based on this approach distinct glycation sites located in vicinity of the LDL-R binding region of ApoB-100 were found in gLDL. Our finding may explain some of the biological dysfunction related to elevated plasma levels of glycosylated LDL species (e.g. in diabetic patients).




Fig. 2 Results of the MALDI Lipid Profiling of PLs isolated from different modified LDL samples. Arrows indicate increase of LPCs (red) and decrease of PUFA-containing PCs (green) as the result of ROS-dependent degradation processes [5].

The study of aging processes using *in vitro* glycation experiments applied to low density lipoproteins using MALDI-QIT-TOF-tandem Mass Spectrometry

Table 3 Component area reduction after 24Hrs treatment of oil in rat plasma




Fig. 3 Positive ion MALDI-TOF-MS/MS spectra of ApoB-100 peptides indicating glycation sites at Lysine 2697 (left) and position 4040 (right). 3 weeks of oxLDL incubation by glucose was needed to detect the modification at AA 4040.




Fig. 4 PTM Finder results based on LC-MALDI-TOF-MS/MS spectra of goxLDL samples searched against a reduced protein database including ApoB-100 and human serum albumin. Known neutral losses related to sugar moieties attached to the AA backbone [7, 8] were found (red circles).

Conclusion

- Our MALDI-QIT-TOF-MS/MS based approach represents a powerful tool to study and monitor the modification process of proteins (e.g. due to the influence of ROS and/or sugars) on the constituent lipid and AA level.
- Software-assisted data analysis using the PTM Finder tool (Shimadzu, UK) allows site-specific tracking of protein modifications (e.g. glycation) within huge MALDI-MS/MS spectra datasets in a very convenient way.
- This approach shows the potential to serve as platform for monitoring of “aging processes”. Such method could easily be translated in food products, pharmaceutical formulations as well as for clinical monitoring (e.g. in diabetic patients) in the future.

Literature

(1) Kaanane *Prog. Clin. Biol. Res.* 1989;304:301-27. (2) Humpert et al. *J. Leuko. Biol.* 2009, 86, 589-597. (3) Witztum and Steinberg *Trends Cardiovasc. Med.* 2001; 11:93-102. (4) Sobal et al. *Prostagl. Leuk. Essent. Fatty Acids* 2000; 63: 177-86. (5) Stübiger et al. *Chem. Phys. Lipids* 2011, 164, 563-572. (6) Stubiger and Belgacem *Anal. Chem.* 2007; 79: 3206-3213. (7) Montgomery et al. *Rapid Commun. Mass Spectrom.* 2010; 24: 1118–1130. (8) Brancia et al. *J. Mass Spectrom.* 2006; 41: 1179–1185.



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.
The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.