

A Semi-Automated Method for Sequencing Oligonucleotides using ISD and Pseudo-MS³ on a MALDI-Ion Trap-TOF Mass Spectrometer

IMSC 2012 PWe-147

Matthew E. Openshaw¹, Omar Belgacem¹ and Marco Smith²

¹Shimadzu, Manchester, UK;

²GlaxoSmithKline, Stevenage, UK

A Semi-Automated Method for Sequencing Oligonucleotides using ISD and Pseudo-MS³ on a MALDI-Ion Trap-TOF Mass Spectrometer

Introduction

Quality Control (QC) analysis of oligonucleotides is commonly performed using MALDI-TOF. Using linear mode analysis, confirmation of the expected mass can be achieved with modest mass accuracy (typically 1-2 Da). An advantage of using MALDI-TOF technology for the analysis of oligonucleotides is the ability to perform in-source decay (ISD) sequencing, if required.

For *unmodified* RNA oligonucleotides, the mass accuracy achieved during ISD analyses is often sufficient to confirm the sequence. However, oligonucleotides designed for use as therapeutics often contain modifications to improve stability and resistance to degradation. In the case of 2'-O-methyl phosphorothioate-modified RNA

oligonucleotides, the mass difference between residues can be as small as 1 Da (e.g. modified C = 336 Da, modified U = 335 Da). In such cases, the lower mass accuracy of linear mode ISD may not be sufficient for unambiguous sequence determination.

Here, we describe an approach using ISD performed on a MALDI-Ion Trap-TOF mass spectrometer. The configuration of this instrument is such that high mass accuracy and monoisotopic resolution are achieved for ISD fragments. In a further development of this application, we applied software originally developed for copolymer analysis, for the semi-automated sequencing of the modified oligonucleotides using MALDI-ISD data.

Experimental

Samples were prepared in deionised water and were desalted (x2) using Dowex 50WX8-200 ion-exchange resin beads (Sigma). Samples were prepared using 3-hydroxypicolinic acid (HPA) (Fluka) and ammonium citrate (Fluka). Sample solution/Dowex resin/HPA matrix/ammonium citrate were mixed in a microcentrifuge tube and the resin beads allowed to settle. An aliquot of this solution was deposited onto a stainless steel MALDI target and dried in a vacuum drier box.

Samples were analysed on an AXIMA Resonance MALDI-Ion Trap-TOF mass spectrometer (Shimadzu, UK). For ISD experiments, the laser power was increased by approx. 10% compared with that used for MS. Samples were analysed in Mid 850 mode (approx. trapped mass range = ~850 – 3500 *m/z*). 500-800 profiles were acquired (2 shots/profile).

MALDI-ISD data were interpreted in a semi-automated manner using *Polymer Analysis* software (Shimadzu). A tolerance of 200 mDa was used when matching candidate oligonucleotide compositions to the experimental data.

Results

Three 2'-O-methyl phosphorothioate-modified (2'-OMe) oligonucleotides samples (labelled 1, 2 and 3) were used to test the proposed workflow. The structures of unmodified and 2'-OMe-modified oligonucleotides differ from

unmodified oligonucleotides in that: (i) one of the non-bridging oxygens in the backbone phosphate is replaced by sulphur and; (ii) the -OH in the 2' position of the nucleoside is replaced with -OCH₃.





Fig. 1 (left) MALDI-Ion Trap-TOF-ISD spectrum obtained for sample 1 (5 mg/mL) and; (right) expanded views showing resolution of selected ISD fragments

A Semi-Automated Method for Sequencing Oligonucleotides using ISD and Pseudo-MS³ on a MALDI-Ion Trap-TOF Mass Spectrometer




Fig. 4 (right): Polymer Analysis software results obtained during semi-automated oligonucleotide sequencing. The upper trace shows the theoretical isotope pattern for the selected composition and the lower trace shows the experimental data.

In Fig. 3 (Polymer Analysis software screenshot), from the iteration shown, 2 compositions are proposed which are valid based on the ratio of 'number of residues' and 'number of phosphorothioate residues' (# residues = # phos groups). Selecting the proposed compositions allows comparison of the theoretical isotope pattern with the experimental data obtained (see Fig. 4). Using the *Polymer Analysis* score (based on isotopic pattern fit), the incorrect composition (2033.2 *m/z* (Fig. 4 (a))) can be quickly eliminated.

The full oligonucleotide sequence was determined by combining sequence information obtained from both the 5'- and 3'-termini. For sample 2, it was not possible to confirm the central 2 residues (shown as XX). Using the intact molecular weight for the sample (determined by

linear mode MALDI-MS, not shown) and the determined 5'- and 3'-sequences, the missing residues were consistent with either CC, [CU] or UU. A higher mass accuracy measurement of the precursor would help eliminate such ambiguities.

For all the samples, the composition of the terminal sequences were proposed in the first iteration but the sequence order was not known. The termini sequences (typically approx. 4 residues) can be confirmed by performing pseudo-MS³ (i.e. MS/MS of the ISD fragment). See Fig. 5). The oligonucleotide samples used in this work were synthetic products and the theoretical sequence was used to facilitate interpretation of the MS/MS spectra for the terminal residues.




Fig. 5: Pseudo-MS³ (i.e. MS/MS of ISD fragments) of terminal fragments of sample 1:
(a) 5'-terminal fragment (1408.1 *m/z*; [UCAA...]) and;
(b) 3'-terminal fragment (1362.1 *m/z*; [...UUCU])

A Semi-Automated Method for Sequencing Oligonucleotides using ISD and Pseudo-MS³ on a MALDI-Ion Trap-TOF Mass Spectrometer

Table 1 shows the sequences determined for samples 1, 2 and 3. The termini sequences (underlined in table 1) were confirmed using the theoretical oligonucleotide sequence and the MS/MS data.

Sample	Sequence obtained	# Correct residues (% correct)
1	<u>UCAAGGAAGAUGGCAUUUCU</u>	20/20 (100 %)
2	<u>GGCCAAACCXXGGCUUCCA</u>	16/20 (80 %)
3	<u>UCAAGGAAGAUGGCAUUUCU</u>	21/21 (100 %)

Table 1: Summary of sequencing results obtained for oligonucleotide samples 1, 2 and 3. For sample 2, the residues shown as XX could not be determined as the higher mass ISD fragment ions were too low abundance.

Note: sequence confirmation of terminal residues (underlined) was performed using the theoretical oligonucleotide sequence.

Conclusion

- The advantages of the MALDI-ion trap-TOF compared with a regular linear mode MALDI-TOF for in-source decay include the high mass accuracy and monoisotopic resolution obtained for the ISD fragments and the ability to perform high quality pseudo-MS³ for terminal sequence confirmation.
- However, the higher mass ISD fragments towards the middle of the sequence have lower abundance which limits the size of the oligonucleotide which can be fully sequenced (approx. 20-mer). Full sequence coverage can be determined by combining sequences from both termini. In the case of sample 2, portions of the sequence were determined incorrectly.
- The *Polymer Analysis* software was shown to be applicable for the semi-automated sequencing of oligonucleotides using MALDI-ion trap-TOF ISD data, although further optimisation is still required. As shown, this software can be used in an iterative approach to determine the sequence of one residue at a time. However, manual verification of the correct residues is still required although this could be automated in the future using simple rules to exclude sequences which are not possible.