


LC/MS, GC/MS Data Analysis Software

# Multi-omics Analysis Package



The Multi-omics Analysis Package, developed for metabolic engineering applications, provides the ability to automatically generate metabolic maps and perform a variety of data analysis for the vast data generated in fields like metabolomics, proteomics and flux analysis. It offers a powerful platform to support drug discovery, bioengineering and other life science research applications.

## Automatically Visualize Quantitative Changes in Metabolites and Proteins on a Metabolic Map

Dramatically reduces the amount of work required for bottleneck processes of analyzing and visualizing data, such as displaying data on a metabolic map and analyzing correlations.


## Works Perfectly with Software Developed by the World's Top Research Institutions

In combination with various other software (called gadgets) connected on the GARUDA™ platform, it offers optimal workflows for multi-dimensional analyses.

## Integrated Support for Operations from Data Acquisition to Analysis

The software includes various data processing gadgets which connect with corresponding data analysis gadgets, providing a seamless, connected experience of a single software program.

It is designed to work in combination with various Shimadzu databases prepared with "Ready-to-Use Methods" for everything from sample pretreatment to analytical conditions, which ensures the entire process flow, from obtaining mass spectrometer measurements to data analysis, can be performed smoothly.



# Data Analysis Example

Changes in metabolite levels in a cell culture medium over time were measured using a GC-MS system. Acquired data was analyzed using the Multi-omics Analysis Package and metabolite changes were visualized by displaying them on a metabolic map.

## Experiment Protocol



Culture MCF-7 cell line for 15 hours.

Replace culture medium with medium spiked with rotenone.

Collect culture medium after 0, 3, 6, 9, and 24 hours.

Use GC-MS system to measure metabolite quantities.

Analyze data using Multi-omics Analysis Package.



## Tools for Data Analysis

The Multi-omics Analysis Package is based on software tools (called gadgets) that have been released on the GARUDA platform - an open research platform, developed by the GARUDA Alliance led by The Systems Biology Institute, Japan (SBI).



<http://www.garuda-alliance.org/>

### Data Analysis Tools Used in the Multi-omics Analysis Package



#### ■ Volcano Plot

A tool that combines a t-test (statistically significant difference) and a fold-change (Example: Difference in mean value such as 2 times or 1/2) to visualize the differences between the two groups. The Volcano Plot gadget developed by Shimadzu is included in the package.



#### ■ VANTED

Tool maintained at University of Konstanz, Germany, for visualization and analysis of networks across different data sets. (GARUDA support was developed at Monash University)



#### ■ iPath

Data analysis tool developed by the European Molecular Biology Laboratory that can be used for visualization of diverse metabolic pathways or data mapping and customization.



#### ■ Cytoscape

Bioinformatics tool developed by the Cytoscape Consortium, used to visualize metabolic pathways, to integrate gene expression profiles with related data, and so on. It is especially useful for analyzing networks and visualizing correlations.

## Method Packages for Omics



### ■ GC/MS, GC-MS/MS Method Package for Metabolites

Description	Flyer code
Smart Metabolites Database™	C146-E277A

### ■ LC/MS/MS Method Packages

Description	Flyer code
Primary Metabolites	C146-E227A
Short Chain Fatty Acid	C146-E355
Lipid Mediators	C146-E381
D/L Amino Acids	C146-E336
Cell Culture Profiling	C146-E279

### ■ LC/MS/MS MRM Library

Description	Flyer code
Metabolic Enzymes in Yeast	C146-E275A
Phospholipid Profiling	C146-E314

Smart Metabolites Database is a trademark of Shimadzu Corporation.  
GARUDA is a trademark of The Systems Biology Institute.



**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.