

Determination of multiple pesticide residues in animal foods using on-line gel permeation chromatography-gas chromatography/mass spectrometry

Pittcon 2014 1680-8P

Ying Ye, Xiaohua Liu, Song Zhan, Taohong Huang,
Shin-ichi Kawano, Yuki Hashi,
Shimadzu Global COE for Application &
Technical Development, Guang Zhou, 510010, China

Determination of multiple pesticide residues in animal foods using on-line gel permeation chromatography-gas chromatography/mass spectrometry

Introduction

An analytical method for the determination of multiple pesticide residues in animal foods by modified QuEChERS method with on-line gel permeation chromatography/gas chromatography/mass spectrometry (GPC/GC/MS). Targets were 21 pesticides: methamidophos, dichlorvos, α -BHC, β -BHC, γ -BHC, δ -BHC, dimethoate, terbufos, chlorpyrifos-methyl, metolachlor, isocarbophos, p, p'-DDE, p, p'-DDD, o, p'-DDT, p, p'-DDT, phentriazophos, permethrin, cypermethrin, flucythrinate, fenvalerate, and deltamethrin. The samples were extracted from homogenized foods with acetonitrile assisted by n-hexane, and separated with liquid-liquid partition. The supernatant liquid was purified by PSA and C18 to remove most of the

fats and pigments in samples, then after on-line GPC-GC/MS analysis which further removed macromolecular interference material, such as fat, the background interference brought about by the complex matrix in samples was effectively reduced. At the spiked level of 0.02 mg/kg, recoveries for most of pesticides were from 72.9% to 118.6%, and the relative standard deviations ranged from 0.88% to 9.67%. The limits of detection and the limits of quantification were 0.6~2.5 μ g/kg and 2.0~8.4 μ g/kg, respectively. The method is simple, rapid and characterized with acceptable sensitivity and accuracy to meet the requirements for the analysis of multiple pesticide residues in animal food.

Experimental

Sample pretreatment




Figure 1 Schematic flow diagram of the sample preparation

Determination of multiple pesticide residues in animal foods using on-line gel permeation chromatography-gas chromatography/mass spectrometry

Instrument

Analytical Conditions

| GPC | |
|------------------|---------------------------------------|
| Mobile phase | : Acetone/Cyclohexane (3/7 by volume) |
| Flow rate | : 0.1mL/min |
| Column | : Shodex CLNpak EV-200 (2mmx150mm) |
| Oven temperature | : 40 °C |
| Injection Volume | : 10 µL |

| GCMS | |
|--------------------------|--|
| Column | : deactivated silica tubing [0.53mm (ID)x5m (L)] +pre-column Rxi-5ms [0.25mm (ID)x5m (L)] Rxi-5ms [0.25mm (ID)x30m (L),Thickness 0.25µm] |
| Injector | : PTV |
| Injector time program | : 120 °C (5min)-(100 °C/min)-250 °C (33.7min) |
| Oven temperature program | : 82 °C (5min) (8 °C/min)-300 °C (8min) |
| Linear velocity | : 45.0cm/sec |
| Ion Source temperature | : 230 °C |
| Interface temperature | : 300 °C |



Figure 2 Shimadzu GPC-GCMS

Determination of multiple pesticide residues in animal foods using on-line gel permeation chromatography-gas chromatography/mass spectrometry

Results




Figure 3 Chromatograms of SIM from mix standards

Table 1 Retention times, linear equation, correlation coefficients, limitis of detection (LODs, $S/N \geq 3$), limitis of quantification (LOQs, $S/N \geq 10$), the averager recoveries and the relative standard deviations (RSDs, $n=3$) for pesticides

| No. | pesticide | t_r (min) | Correlation Coefficient* | LOD ($\mu\text{g}/\text{kg}$) | LOQ ($\mu\text{g}/\text{kg}$) | pork | | shrimp | |
|-----|---------------------|----------------|-----------------------------|------------------------------------|------------------------------------|--------------|---------|--------------|---------|
| | | | | | | Recovery (%) | RSD (%) | Recovery (%) | RSD (%) |
| 1 | Methamidophos | 11.800 | 0.9993 | 1.8 | 6.0 | 72.9 | 2.99 | 83.4 | 4.52 |
| 2 | Dichlorvos | 11.983 | 0.9999 | 0.6 | 2.0 | 106.3 | 4.47 | 106.6 | 2.87 |
| 3 | α -BHC | 19.742 | 0.9999 | 0.8 | 2.7 | 97.3 | 8.22 | 91.4 | 6.83 |
| 4 | Dimethoatee | 20.108 | 0.9995 | 1.8 | 6.0 | 98.4 | 7.79 | 96.3 | 6.87 |
| 5 | β -BHC | 20.475 | 0.9999 | 1.7 | 5.7 | 95.6 | 6.80 | 97.8 | 6.43 |
| 6 | γ -BHC | 20.708 | 0.9999 | 0.8 | 2.7 | 105.2 | 4.22 | 98.0 | 6.15 |
| 7 | Terbufos | 20.817 | 0.9999 | 0.8 | 2.7 | 99.3 | 3.52 | 95.5 | 6.42 |
| 8 | δ -BHC | 21.433 | 0.9999 | 1.7 | 5.7 | 119.4 | 9.67 | 100.0 | 6.70 |
| 9 | Chlorpyrifos-methyl | 22.358 | 0.9999 | 1.6 | 5.4 | 118.6 | 9.52 | 98.3 | 3.58 |
| 10 | Metolachlor | 23.567 | 0.9999 | 1.6 | 5.4 | 110.4 | 3.71 | 95.8 | 4.71 |
| 11 | Isocarbopho | 23.883 | 0.9993 | 2.1 | 7.0 | 110.7 | 3.21 | 93.1 | 5.11 |
| 12 | p, p'-DDE | 26.192 | 0.9999 | 0.8 | 2.7 | 73.0 | 0.88 | 70.1 | 4.98 |
| 13 | p, p'-DDD | 27.250 | 0.9999 | 0.8 | 2.7 | 106.4 | 2.67 | 103.4 | 3.62 |
| 14 | o, p'-DDT | 27.342 | 0.9998 | 0.8 | 2.7 | 86.5 | 7.88 | 87.9 | 7.25 |
| 15 | Phentriazophos | 27.575 | 0.9996 | 1.0 | 3.4 | 108.6 | 4.25 | 98.7 | 6.07 |
| 16 | p, p'-DDT | 28.192 | 0.9998 | 1.6 | 5.4 | 99.0 | 1.32 | 100.9 | 6.93 |
| 17 | Permethrin-1 | 31.617 | 0.9998 | 1.8 | 6.0 | 93.4 | 3.85 | 90.8 | 4.98 |
| 18 | Permethrin-2 | 31.783 | | | | | | | |
| 19 | Cypermethrin-1 | 32.767 | 0.9996 | 1.9 | 6.4 | 98.4 | 3.85 | 95.4 | 5.80 |
| 20 | Cypermethrin-2 | 32.900 | | | | | | | |
| 21 | Cypermethrin-3 | 33.000 | | | | | | | |
| 22 | Cypermethrin-4 | 33.050 | | | | | | | |
| 23 | Flucythrinate-1 | 33.000 | 0.9996 | 1.8 | 6.0 | 102.6 | 1.97 | 88.6 | 6.25 |
| 24 | Flucythrinate-2 | 33.275 | | | | | | | |
| 25 | Fenvalerate-1 | 34.167 | 0.9997 | 1.8 | 6.0 | 104.7 | 3.50 | 101.1 | 7.25 |
| 26 | Fenvalerate-1 | 34.500 | | | | | | | |
| 27 | Deltamethrin | 33.475 | 0.9988 | 2.5 | 8.4 | 106.1 | 3.42 | 93.3 | 5.75 |

* Concentration range: 10 $\mu\text{g}/\text{L}$ ~500 $\mu\text{g}/\text{L}$

Determination of multiple pesticide residues in animal foods using on-line gel permeation chromatography-gas chromatography/mass spectrometry

Conclusion

A based on automated online gel permeation chromatography coupled with gas chromatography/mass spectrometry (GPC-GC/MS) method was developed for the determination of multiple pesticide residues in animal

foods by QuEchERS method. So provides a simple, rapid, and reliable method for the detection of pesticide residues in animal foods.