

Shimadzu LC Column Kits for Shimadzu HPLC/UHPLC Systems

Shim-pack Method Development Column Kit

CoreFocus



LC column kits for optimal column selection

Kit types	HPLC	UHPLC	HPLC (LC-MS)	UHPLC (LC-MS)
① L1 Kit C18 only	0			
② L1 Kit for HPLC / UHPLC (LC-MS)	0	0	0	0
③ Maximum Selectivity RP Kit for HPLC / UHPLC Type A	0	0	0	0
Maximum Selectivity RP Kit for HPLC / UHPLC Type B	0	0	0	0
Maximum Selectivity RP Kit for HPLC / UHPLC (LC-MS)	0	0	0	0

^{*}These column kits do not guarantee the appropriate separation for customer analyses.

○: Most suitable

O: Compatible

■Kit composition columns

① L1 Kit for HPLC						
For HPLC (150 mm × 4.6 mm	For HPLC (150 mm × 4.6 mm I.D., 5 μm)					
Scepter C18-120	227-31020-05					
Scepter HD-C18-80	227-31024-05					
Arata C18	227-32805-04	227-32810-01				
GISS C18	227-30061-06	227-32010-01				
GIST C18-AQ	227-30742-07					
Velox C18	227-32012-03					

② L1 Kit for HPLC / UHPLC (LC-MS)					
For HPLC (150 mm × 4.6 mm l.D., 5 μm)		For UHPLC (100 mm × 2.1 mm I.D., 1.8 / 1.9 / 2 μm)			
Scepter C18-120	227-31020-05	227-32810-02	Scepter C18-120	227-31012-05	
GISS C18	227-30061-06		GISS C18	227-30048-02	
GIST C18	227-30017-07		GIST C18	227-30001-04	227-32810-03
GIST C18-AQ	227-30742-07		GIST C18-AQ	227-30807-02	227-32610-03
Velox C18	227-32012-03		Velox C18	227-32007-03	
Velox SP-C18	227-32006-03		Velox SP-C18	227-32001-03	

③ Maximum Selectivity RP Kit for HPLC / UHPLC Type A					
For HPLC (150 mm × 4.6 mm I.D., 5 μm)		For UHPLC (100 mm × 2.1 mm I.D., 1.8 / 1.9 μm)			
Scepter C18-120	227-31020-05	- 227-32810-04	Scepter C18-120	227-31012-05	
Scepter HD-C18-80	227-31024-05		Scepter HD-C18-80	227-31026-05	
Scepter C8-120	227-31041-05		Scepter C8-120	227-31033-05	227-32810-05
Scepter Phenyl-120	227-31071-05		Scepter Phenyl-120	227-31063-05	227-32810-03
Scepter PFPP-120	227-31061-05		Scepter PFPP-120	227-31053-05	
Velox Biphenyl	227-32018-03		Velox Biphenyl	227-32013-03	


Maximum Selectivity RP Kit for HPLC / UHPLC Type B					
For HPLC (150 mm × 4.6 mm I.D., 5 μm)		For UHPLC (100 mm × 2.1 mm I.D., 1.8 / 1.9 / 2 / 3 μm)			
GIST C18	227-30017-07	- 227-32810-06	GIST C18	227-30001-04	
GIST C18-AQ	227-30742-07		GIST C18-AQ	227-30807-02	
GIST C8	227-30173-07		GIST C8	227-30160-04	227-32810-07
GIST Phenyl-Hexyl	227-30690-05		GIST Phenyl-Hexyl	227-30713-03	227-32610-07
Velox PFPP	227-32024-03		Velox PFPP	227-32019-03	
Velox Biphenyl	227-32018-03		Velox Biphenyl	227-32013-03	

⑤ Maximum Selectivity RP Kit for HPLC / UHPLC (LC-MS)					
For HPLC (150 mm × 4.6 mm I.D., 5 μm)		For UHPLC (100 mm × 2.1 mm I.D., 1.8 / 1.9 μm)			
Scepter C18-120	227-31020-05	- 227-32810-08	Scepter C18-120	227-31012-05	
GIST C18-AQ	227-30742-07		GIST C18-AQ	227-30807-02]
GISS C18	227-30061-06		GISS C18	227-30048-02	227-32810-09
Velox C18	227-32012-03		Velox C18	227-32007-03	227-32810-03
Velox PFPP	227-32024-03		Velox PFPP	227-32019-03	1
Velox Biphenyl	227-32018-03		Velox Biphenyl	227-32013-03	


1 L1 Kit For HPLC

Peaks 1. Uracil 5. o-Terphenyl 2. Caffeine 6. Amylbenzene 3. Phenol 7. Triphenylene 4. Butylbenzene


Shim-pack Scepter™ C18-120




Shim-pack Arata™ C18




Shim-pack Scepter HD-C18-80




Shim-pack™ GISS C18



Shim-pack GIST C18-AQ







② L1 Kit for HPLC / UHPLC (LC-MS)

1. Uracil 5. o-Terphenyl
2. Caffeine 6. Amylbenzene
3. Phenol 7. Triphenylene
4. Butylbenzene

Shim-pack Scepter C18-120




Shim-pack GISS C18




reversed-phase nucleic acid drug analysis Shim-pack GIST C18

*The P/N in this brochure has a stainless body.


An inert body (Claris) is the first choice in




Shim-pack Velox C18



Shim-pack GIST C18-AQ




Shim-pack Velox SP-C18




3 Maximum Selectivity RP Kit for HPLC / UHPLC Type A

Peaks 1. Uracil 5. o-Terphenyl 2. Caffeine 6. Amylbenzene 3. Phenol 7. Triphenylene 4. Butylbenzene


Shim-pack Scepter C18-120




Shim-pack Scepter Phenyl-120




Shim-pack Scepter HD-C18-80




Shim-pack Scepter PFPP-120



Shim-pack Scepter C8-120



Shim-pack Velox Biphenyl




4 Maximum Selectivity RP Kit for HPLC / UHPLC Type B


Peaks

1. Uracil 5. o-Terphenyl
2. Caffeine 6. Amylbenzene
3. Phenol 7. Triphenylene
4. Butylbenzene


Shim-pack GIST C18




Shim-pack GIST Phenyl-Hexyl




Shim-pack GIST C18-AQ




Shim-pack Velox PFPP



Shim-pack GIST C8







5 Maximum Selectivity RP Kit for HPLC / UHPLC (LC-MS)

Peaks 1. Uracil 5. o-Terphenyl 2. Caffeine 6. Amylbenzene 3. Phenol 7. Triphenylene 4. Butylbenzene


Shim-pack Scepter C18-120




Shim-pack Velox C18




Shim-pack GIST C18-AQ




Shim-pack Velox PFPP



Shim-pack GISS C18



Shim-pack Velox Biphenyl



CoreFocus, Shim-pack, Shim-pack Arata, Shim-pack Scepter, and Shim-pack Velox are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



Shimadzu Corporation www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.